



SLC34A3 Antibody, Biotin conjugated

Product Code	CSB-PA854015LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q8N130
Immunogen	Recombinant Human Sodium-dependent phosphate transport protein 2C protein (1-76AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium-dependent phosphate transport protein 2C (Sodium-phosphate transport protein 2C) (Na(+)-dependent phosphate cotransporter 2C) (Sodium/inorganic phosphate cotransporter IIC) (Sodium/phosphate cotransporter 2C) (Na(+)/Pi cotransporter 2C) (NaPi-2c) (Solute carrier family 34 member 3), SLC34A3, NPT2C NPTIIC
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	SLC34A3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.