



MICU1 Antibody, FITC conjugated

Product Code	CSB-PA860652LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BPX6
Immunogen	Recombinant Human Calcium uptake protein 1, mitochondrial protein (55-191AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Calcium uptake protein 1, mitochondrial (Atopy-related autoantigen CALC) (ara CALC) (Calcium-binding atopy-related autoantigen 1) (allergen Hom s 4), MICU1, CALC CBARA1
Immunogen Species	Homo sapiens (Human)
Research Area	Immunology
Target Names	MICU1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.