



# PIN4 Antibody, Biotin conjugated

|                            |                                                                                                                                                                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>        | CSB-PA861478LD01HU                                                                                                                                                                                                                                                  |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                                                                                                                                                                                                       |
| <b>Uniprot No.</b>         | Q9Y237                                                                                                                                                                                                                                                              |
| <b>Immunogen</b>           | Recombinant Human Peptidyl-prolyl cis-trans isomerase NIMA-interacting 4 protein (1-131AA)                                                                                                                                                                          |
| <b>Raised In</b>           | Rabbit                                                                                                                                                                                                                                                              |
| <b>Species Reactivity</b>  | Human                                                                                                                                                                                                                                                               |
| <b>Tested Applications</b> | ELISA                                                                                                                                                                                                                                                               |
| <b>Form</b>                | Liquid                                                                                                                                                                                                                                                              |
| <b>Conjugate</b>           | Biotin                                                                                                                                                                                                                                                              |
| <b>Storage Buffer</b>      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4                                                                                                                                                                                    |
| <b>Purification Method</b> | >95%, Protein G purified                                                                                                                                                                                                                                            |
| <b>Isotype</b>             | IgG                                                                                                                                                                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                                                                                                                                                                          |
| <b>Alias</b>               | Peptidyl-prolyl cis-trans isomerase NIMA-interacting 4 (EC 5.2.1.8) (Parvulin-14) (Par14) (hPar14) (Parvulin-17) (Par17) (hPar17) (Peptidyl-prolyl cis-trans isomerase Pin4) (PPlase Pin4) (Peptidyl-prolyl cis/trans isomerase EPVH) (hEPVH) (Rotamase Pin4), PIN4 |
| <b>Immunogen Species</b>   | Homo sapiens (Human)                                                                                                                                                                                                                                                |
| <b>Research Area</b>       | Epigenetics and Nuclear Signaling                                                                                                                                                                                                                                   |
| <b>Target Names</b>        | PIN4                                                                                                                                                                                                                                                                |
| <b>Usage</b>               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.                                                                                                                                                                                         |