



TRPV6 Antibody, Biotin conjugated

Product Code	CSB-PA861992LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1D0
Immunogen	Recombinant Human Transient receptor potential cation channel subfamily V member 6 protein (620-764AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Transient receptor potential cation channel subfamily V member 6 (TrpV6) (CaT-like) (CaT-L) (Calcium transport protein 1) (CaT1) (Epithelial calcium channel 2) (ECaC2), TRPV6, ECAC2
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	TRPV6
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.