



PCGF1 Antibody, Biotin conjugated

Product Code	CSB-PA863114LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSM1
Immunogen	Recombinant Human Polycomb group RING finger protein 1 protein (113-187AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycomb group RING finger protein 1 (Nervous system Polycomb-1) (NSPc1) (RING finger protein 68), PCGF1, NSPC1 RNF68
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PCGF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.