



PCGF1 Antibody, Biotin conjugated

Product Code	CSB-PA863114LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BSM1
Immunogen	Recombinant Human Polycomb group RING finger protein 1 protein (113-187AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Polycomb group RING finger protein 1 (Nervous system Polycomb-1) (NSPc1) (RING finger protein 68), PCGF1, NSPC1 RNF68
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	PCGF1