



RAB11FIP5 Antibody, Biotin conjugated

Product Code	CSB-PA863141LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BXF6
Immunogen	Recombinant Human Rab11 family-interacting protein 5 protein (374-545AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Rab11 family-interacting protein 5 (Rab11-FIP5) (Gamma-SNAP-associated factor 1) (Gaf-1) (Phosphoprotein pp75) (Rab11-interacting protein Rip11), RAB11FIP5, GAF1 KIAA0857 RIP11
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	RAB11FIP5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.