



RASSF7 Antibody, Biotin conjugated

Product Code	CSB-PA863596LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q02833
Immunogen	Recombinant Human Ras association domain-containing protein 7 protein (300-368AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ras association domain-containing protein 7 (HRAS1-related cluster protein 1), RASSF7, C11orf13 HRC1
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	RASSF7
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.