



WNT8A Antibody, Biotin conjugated

Product Code	CSB-PA863951LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H1J5
Immunogen	Recombinant Human Protein Wnt-8a protein (223-309AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protein Wnt-8a (Protein Wnt-8d), WNT8A, WNT8D
Immunogen Species	Homo sapiens (Human)
Research Area	Stem Cells
Target Names	WNT8A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.