



LANCL2 Antibody, Biotin conjugated

Product Code	CSB-PA865129LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NS86
Immunogen	Recombinant Human LanC-like protein 2 protein (2-106AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	LanC-like protein 2 (Testis-specific adriamycin sensitivity protein), LANCL2, GPR69B TASP
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	LANCL2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.