



PCDH11X Antibody, Biotin conjugated

Product Code	CSB-PA866326LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BZA7
Immunogen	Recombinant Human Protocadherin-11 X-linked protein (987-1117AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Protocadherin-11 X-linked (Protocadherin-11) (Protocadherin on the X chromosome) (PCDH-X) (Protocadherin-S), PCDH11X, KIAA1326 PCDH11 PCDHX
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	PCDH11X
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.