



PIGV Antibody, HRP conjugated

Product Code	CSB-PA868321LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NUD9
Immunogen	Recombinant Human GPI mannosyltransferase 2 protein (400-469AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	GPI mannosyltransferase 2 (EC 2.4.1.-) (GPI mannosyltransferase II) (GPI-MT-II) (Phosphatidylinositol-glycan biosynthesis class V protein) (PIG-V), PIGV
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	PIGV
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.