



KLK14 Antibody, Biotin conjugated

Product Code	CSB-PA868375LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0G3
Immunogen	Recombinant Human Kallikrein-14 protein (114-225AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kallikrein-14 (hK14) (EC 3.4.21.-) (Kallikrein-like protein 6) (KLK-L6), KLK14, KLKL6
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	KLK14
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.