



CEP250 Antibody, FITC conjugated

Product Code	CSB-PA871581LC01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BV73
Immunogen	Recombinant Human Centrosome-associated protein CEP250 protein (40-293AA)
Raised In	Rabbit
Species Reactivity	Human
Form	Liquid
Conjugate	FITC
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Centrosome-associated protein CEP250 (250 kDa centrosomal protein) (Cep250) (Centrosomal Nek2-associated protein 1) (C-Nap1) (Centrosomal protein 2), CEP250, CEP2 CNAP1
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	CEP250
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.