



KIFC1 Antibody, Biotin conjugated

Product Code	CSB-PA871587LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BW19
Immunogen	Recombinant Human Kinesin-like protein KIFC1 protein (337-498AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Kinesin-like protein KIFC1 (Kinesin-like protein 2) (Kinesin-related protein HSET), KIFC1, HSET KNSL2
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KIFC1