



KCMF1 Antibody, HRP conjugated

Product Code	CSB-PA873718LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0J7
Immunogen	Recombinant Human E3 ubiquitin-protein ligase KCMF1 protein (98-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase KCMF1 (EC 2.3.2.27) (FGF-induced in gastric cancer) (Potassium channel modulatory factor) (PCMF) (RING-type E3 ubiquitin transferase KCMF1) (ZZ-type zinc finger-containing protein 1), KCMF1, FIGC ZZZ1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KCMF1