



KCMF1 Antibody, FITC conjugated

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA873718LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9P0J7 |
| Immunogen | Recombinant Human E3 ubiquitin-protein ligase KCMF1 protein (98-226AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | E3 ubiquitin-protein ligase KCMF1 (EC 2.3.2.27) (FGF-induced in gastric cancer) (Potassium channel modulatory factor) (PCMF) (RING-type E3 ubiquitin transferase KCMF1) (ZZ-type zinc finger-containing protein 1), KCMF1, FIGC ZZZ1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Signal Transduction |
| Target Names | KCMF1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |