



KCMF1 Antibody, Biotin conjugated

Product Code	CSB-PA873718LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9P0J7
Immunogen	Recombinant Human E3 ubiquitin-protein ligase KCMF1 protein (98-226AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	E3 ubiquitin-protein ligase KCMF1 (EC 2.3.2.27) (FGF-induced in gastric cancer) (Potassium channel modulatory factor) (PCMF) (RING-type E3 ubiquitin transferase KCMF1) (ZZ-type zinc finger-containing protein 1), KCMF1, FIGC ZZZ1
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	KCMF1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.