



FERMT1 Antibody, FITC conjugated

| | |
|----------------------------|---|
| Product Code | CSB-PA874786LC01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9BQL6 |
| Immunogen | Recombinant Human Fermitin family homolog 1 protein (321-420AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Form | Liquid |
| Conjugate | FITC |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Fermitin family homolog 1 (Kindlerin) (Kindlin syndrome protein) (Kindlin-1) (Unc-112-related protein 1), FERMT1, C20orf42 KIND1 URP1 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Tags & Cell Markers |
| Target Names | FERMT1 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |