



FERMT1 Antibody, Biotin conjugated

Product Code	CSB-PA874786LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9BQL6
Immunogen	Recombinant Human Fermitin family homolog 1 protein (321-420AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Fermitin family homolog 1 (Kindlerin) (Kindlin syndrome protein) (Kindlin-1) (Unc-112-related protein 1), FERMT1, C20orf42 KIND1 URP1
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	FERMT1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.