



MMADHC Antibody, Biotin conjugated

Product Code	CSB-PA875665LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H3L0
Immunogen	Recombinant Human Methylmalonic aciduria and homocystinuria type D protein, mitochondrial protein (26-142AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Methylmalonic aciduria and homocystinuria type D protein, mitochondrial (CbID), MMADHC, C2orf25 CL25022
Immunogen Species	Homo sapiens (Human)
Research Area	Tags & Cell Markers
Target Names	MMADHC
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.