



NAA60 Antibody, Biotin conjugated

Product Code	CSB-PA875691LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9H7X0
Immunogen	Recombinant Human N-alpha-acetyltransferase 60 protein (1-94AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	N-alpha-acetyltransferase 60 (hNaa60) (EC 2.3.1.259) (Histone acetyltransferase type B protein 4) (HAT4) (EC 2.3.1.48) (N-acetyltransferase 15) (N-alpha-acetyltransferase F) (NatF), NAA60, HAT4 NAT15
Immunogen Species	Homo sapiens (Human)
Research Area	Cell Biology
Target Names	NAA60
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.