




AMBN Antibody

Product Code	CSB-PA882061LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NP70
Immunogen	Recombinant Human Ameloblastin protein (264-387AA)
Raised In	Rabbit
Species Reactivity	Human, Rat
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Ameloblastin, AMBN
Immunogen Species	Homo sapiens (Human)
Research Area	Others
Target Names	AMBN

Image



Western Blot
 Positive WB detected in: Rat brain tissue
 All lanes: AMBN antibody at 4µg/ml
 Secondary
 Goat polyclonal to rabbit IgG at 1/50000 dilution
 Predicted band size: 49, 47 kDa
 Observed band size: 49 kDa

Usage For Research Use Only. Not for use in diagnostic or therapeutic procedures.