



SCN11A Antibody, HRP conjugated

Product Code	CSB-PA883417LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UI33
Immunogen	Recombinant Human Sodium channel protein type 11 subunit alpha protein (403-551AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Sodium channel protein type 11 subunit alpha (Peripheral nerve sodium channel 5) (PN5) (Sensory neuron sodium channel 2) (Sodium channel protein type XI subunit alpha) (Voltage-gated sodium channel subunit alpha Nav1.9) (hNaN), SCN11A, SCN12A SNS2
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	SCN11A
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.