



HOXD1 Antibody, HRP conjugated

Product Code	CSB-PA884426LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9GZZ0
Immunogen	Recombinant Human Homeobox protein Hox-D1 protein (59-131AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Homeobox protein Hox-D1 (Homeobox protein Hox-GG), HOXD1, HOX4 HOX4G
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	HOXD1
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.