



MBIP Antibody, Biotin conjugated

Product Code	CSB-PA889101LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9NS73
Immunogen	Recombinant Human MAP3K12-binding inhibitory protein 1 protein (91-241AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	MAP3K12-binding inhibitory protein 1 (MAPK upstream kinase-binding inhibitory protein) (MUK-binding inhibitory protein), MBIP
Immunogen Species	Homo sapiens (Human)
Research Area	Signal Transduction
Target Names	MBIP
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.