



PPIE Antibody, HRP conjugated

| | |
|----------------------------|--------------------------------------------------------------------------------------------------------------------------|
| Product Code | CSB-PA890944LB01HU |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9UNP9 |
| Immunogen | Recombinant Human Peptidyl-prolyl cis-trans isomerase E protein (86-141AA) |
| Raised In | Rabbit |
| Species Reactivity | Human |
| Tested Applications | ELISA |
| Form | Liquid |
| Conjugate | HRP |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | >95%, Protein G purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Alias | Peptidyl-prolyl cis-trans isomerase E (PPIase E) (EC 5.2.1.8) (Cyclophilin E) (Cyclophilin-33) (Rotamase E), PPIE, CYP33 |
| Immunogen Species | Homo sapiens (Human) |
| Research Area | Epigenetics and Nuclear Signaling |
| Target Names | PPIE |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |