



NUDT5 Antibody, Biotin conjugated

Product Code	CSB-PA891535LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9UKK9
Immunogen	Recombinant Human ADP-sugar pyrophosphatase protein (34-166AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	ADP-sugar pyrophosphatase (EC 3.6.1.13) (8-oxo-dGDP phosphatase) (EC 3.6.1.58) (Nuclear ATP-synthesis protein NUDIX5) (EC 2.7.7.96) (Nucleoside diphosphate-linked moiety X motif 5) (Nudix motif 5) (hNUDT5) (YSA1H), NUDT5, NUDIX5
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	NUDT5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.