



ZKSCAN5 Antibody, Biotin conjugated

Product Code	CSB-PA897093LD01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y2L8
Immunogen	Recombinant Human Zinc finger protein with KRAB and SCAN domains 5 protein (490-643AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	Biotin
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Zinc finger protein with KRAB and SCAN domains 5 (Zinc finger protein 95 homolog) (Zfp-95), ZKSCAN5, KIAA1015 ZFP95
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	ZKSCAN5
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.