



NUMBL Antibody, HRP conjugated

Product Code	CSB-PA897595LB01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9Y6R0
Immunogen	Recombinant Human Numb-like protein (447-576AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA
Form	Liquid
Conjugate	HRP
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Numb-like protein (Numb-related protein) (Numb-R), NUMBL
Immunogen Species	Homo sapiens (Human)
Research Area	Neuroscience
Target Names	NUMBL
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.