




EGR3 Antibody

Product Code	CSB-PA007486LA01HU
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q06889
Immunogen	Recombinant Human Early growth response protein 3 protein (201-281AA)
Raised In	Rabbit
Species Reactivity	Human
Tested Applications	ELISA, WB; Recommended dilution: WB:1:500-1:5000
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	>95%, Protein G purified
Isotype	IgG
Clonality	Polyclonal
Alias	Early growth response protein 3 (EGR-3) (Zinc finger protein pilot), EGR3, PILOT
Immunogen Species	Homo sapiens (Human)
Research Area	Epigenetics and Nuclear Signaling
Target Names	EGR3

Image



Western Blot

Positive WB detected in: MCF-7 whole cell lysate, A549 whole cell lysate

All lanes: EGR3 antibody at 3.5µg/ml

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 43, 39 kDa

Observed band size: 43 kDa