



Recombinant Lactobacillus plantarum Glycerol uptake facilitator protein-like 6 (glpF6)

Product Code	CSB-CF950889LMS
Uniprot No.	F9UUD2
Product Type	Transmembrane Protein
Immunogen Species	Lactobacillus plantarum (strain ATCC BAA-793 / NCIMB 8826 / WCFS1)
Sequence	MRKYLAEFLGTFMLVFLGTATVVIKGDVLAIGLAFGLAITVSAYAFGGISGGHF NPAVTTAMLINRRIDAADAIGYIIAQIIGAIVASAAVKSFVSALGLSATLLGQTDFP KIGSGMAFFVEALVTFLFLMVILNVTSNDHGNADFAGLTIGVTLAFLIIVALNLTG GSLNPARSIGPAIFAGGSALSHLWVYILAPEVGAILAAFCARVMGSED
Source	in vitro E.coli expression system
Target Names	glpF6
Expression Region	1-216
Protein Length	Full Length
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.