



Recombinant Lactobacillus plantarum D/L-lactic acid transporter (IarD)

Product Code	CSB-CF167360LMS
Uniprot No.	F9UST3
Product Type	Transmembrane Protein
Immunogen Species	Lactobacillus plantarum (strain ATCC BAA-793 / NCIMB 8826 / WCFS1)
Sequence	MVHQLIAEFMGTALMIIFGVGVHCSSVLKGTKYRGS GHIFAITTWGFGISVALFI FGNVCINPAMVLAQCLLGNIAWSLFIPYSVAEVLGGVVGSVIVWIMYADHFKAS TDEISPITIRNLFCTAPAVRNLPNFFVELFDTFIFISGILAISEIKTPGIVPIGVLL VWAIGMGLGGPTGFAMNLRDMGPRIAHAILPIANKADSDWQYGIIVPGIAPFV GAAIAAWFMHGFFGIN
Source	in vitro E.coli expression system
Target Names	IarD
Expression Region	1-238
Protein Length	Full Length
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.