



Recombinant Rat Alpha-1,6-mannosylglycoprotein 6-beta-N-acetylglucosaminyltransferase A (Mgat5)

Product Code	CSB-CF147087RA
Uniprot No.	Q08834
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	MAFFSPWKLSSQKLGFFLVTFGFIWGMMLLHFTIQQRTQPESSSMLREQILDLSKRYIKALAEENRNVDGYPYAGVMTAYDLKKTAVLLDNILQRIGKLESKVDNLVNGTGANSTNSTTAVPSLVSLEKINVADIINGVQEKCVLPPMDGYPHCEGKIKWMKDMWRS DPCYADYGVDGTSCSFFIYLSEVENWCPRLPWRANKNPYEEADHNSLAEIRTD FNILY GMMKKHEEFRWMRLRIRRMADAWIQAIKSLAEKQNLEKRKRKKILVHLGLLTKESGFKIAETA FSGGPLGELVQWSDLITSLYLLGHDIRISASLAELKEIMKKVVGNRSGCPTVGDRIVELIYIDIVGLAQFKKTLGPSWVHYQCMLRVLDSFGTEPEFNHASYAQSKGHKTPWGKWNLNPPQQFYTMFPHTPDNSFLGFVVEQHLNSSDIHHINEIKRQNSLVYGVKVD SFWKNKKIYLDIIHTYMEVHATVYGSSTKNIPSYVKNHGILSGRDLQFLLRETKLFVGLGFPYEGPAPLEAIANGCAFLNPKFNPPKSSKNTDFFIGKPTLRELTSQHPYAEVFIGRPHVWTVDLNNREEVEDAVKAILNQKIEPYMPYEFTCEGMLQRINAFIEKQDFCHGQVMWPPLSALQVKLAEPGQSKQVCQESQLICEPSFFQHLNKEKDLLKYKVICQSSELYKDILVPSFYPKSKHCVFQGDLLL FSCAGAHPTHQRICPCRDFIKGQVALCKDCL
Source	in vitro E.coli expression system
Target Names	Mgat5
Expression Region	1-740
Protein Length	Full Length
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.