



# Recombinant Rat Ammonium transporter Rh type A (Rhag)

<b>Product Code</b>	CSB-CF756017RA
<b>Uniprot No.</b>	Q7TNK7
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Rattus norvegicus (Rat)
<b>Sequence</b>	MRFKFSLIALSLEVVMIVSFALFVEYETSQNGSQKSASQQNASQQNAAAQQNA SQQGNASSPAKEDQFFQLYPLFQHVHVMIFVGFGLMTFLKKYGFSGVGFNL FLAALGLQWGTIVQGLLHSHGLKFPFRIKNMINADFSTATVLISFGAVLGKTSPI QMIIMTILEIAVFAGNEHLVTEIFKASDTGASMTIHAFGAYFGLAVAGVLYRSGLK HGHPNEESVYHSDLFAMIGTLFLWMFWPSFNSAIAQPENNQYRAIVNTYMSLA ACVITAYALSSLVERRGRLDMVHIQNATLAGGVAVGTCADMEIPLYFAMTIGSIA GIISVLGYKFLSPLLAHKLMIHDTCGVHNLHGLPGVFGGLASIVAIWGWKSTVST MAMQATALGSSIGSAIVGGLVTGLILKLPVWNQPPDEYCFDDSVSWKVPKYRE LDNYFFQHVTNHVEHEV
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	Rhag
<b>Expression Region</b>	1-450
<b>Protein Length</b>	Full Length
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.