



# Recombinant Human Solute carrier family 40 member 1 (SLC40A1)

<b>Product Code</b>	CSB-CF865089HU
<b>Uniprot No.</b>	Q9NP59
<b>Product Type</b>	Transmembrane Protein
<b>Immunogen Species</b>	Homo sapiens (Human)
<b>Sequence</b>	MTRAGDHNRRQGCCGSLADYLTSAKFLLYLGHSLSTWGDRMWHFAVSVFLV ELYGNSLLLTAVYGLVVAGSVLVLGAIIGDWVDKNARLKVAQTSLVVQNVSVIL CGIILMMVFLHKHELLTMYHGWVLTSCYILIITIANIANLASTATAITIQRDWIVVVA GEDRSKLANMNATIRRIDQLTNILAPMAVGQIMTFGSPVIGCGFISGWNLVSMC VEYVLLWKVYQKTPALAVKAGLKEEETELKQLNLHKDTEPKPLEGTHLMGVKD SNIHELEHEQEPTCASQMAEPFRFTFRDGVWSYYPVFLAGMGLAFLYMTVL GFDCITTGYAYTQGLSGSILMGSASAITGIMGTVAFTWLRRKCGLVRTGLISG LAQLSCLILCVISVFMPSPLDLSVSPFEDIRSRFIQGESITPTKIPEITTEIYMSN GSNSANIVPETSPEVPIISVLLFAGVIAARIGLWSFDLTVTQLLQENVIESERG IINGVQNSMNYLLDLLHFIMVILAPNPEAFGLLVLSVSFVAMGHIMYFRFAQNTL GNKLFACGPDAKEVRKENQANTSVV
<b>Source</b>	in vitro E.coli expression system
<b>Target Names</b>	SLC40A1
<b>Expression Region</b>	1-571
<b>Protein Length</b>	Full Length
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.