



Recombinant Rat Neuronal acetylcholine receptor subunit alpha-3 (Chrna3)

Product Code	CSB-CF005389RA
Uniprot No.	P04757
Product Type	Transmembrane Protein
Immunogen Species	Rattus norvegicus (Rat)
Sequence	SEAEHRLFQYLFEDYNEIIRPVANVSHPVIIQFEVSMSQLVKVDEVNQIMETNL WLKQIWNDYKWKWSPDYQGVEFMRVPAEKIWKPDIVLYNNADGDFQVDDKT KALLKYTG EVTWIPPAIFKSSCKIDVTYFPFDYQNMKFGSWSYDKAKIDLVI GSSMNLKDYWESGEWAIKAPGYKHEIKYNCCEEIYQDITYSLYIRRLPLFYTIN LIIPCLLISFLTVLVLYLPDCGKVTLCISVLLSLTVFLLVITETIPSTSLVIPLIGEY LLFTMIFVTL SIVITVFVLNVHYRTPHTHTMPTWVKAVFLNLLPRVMFMTRPTSG EGDTPKTRTFYGAELSNLNCFSRADSKSCKEGYPCQDGTGCGYCHHRRVKISN FSANLTRSSSES VNAVLSLSALSPEIKEAIQSVKYIAENMKAQNVAKEIQDDW KYVAMVIDRIFLWVFLVCILGTAGLFLQPLMARDDT
Source	in vitro E.coli expression system
Target Names	Chrna3
Expression Region	26-499
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.