



Recombinant *Saccharomyces cerevisiae* Aquaporin-like protein 2 (AQY2-2)

Product Code	CSB-CF409777STA
Uniprot No.	A7A0K8
Product Type	Transmembrane Protein
Immunogen Species	<i>Saccharomyces cerevisiae</i> (strain YJM789) (Baker's yeast)
Sequence	MSNESNDLEKNISHLDPTGVDNAYIPPEQPETKHSRNFNIDRDLRNHFIAAVGE FCGTFMFLWCAYVICNVANHDVALTTEPEGSHPGQLIMIALGFGFSVMFSIWC FWWGFEPSRFLFVFGQSHLTSQMCSDDVSSDHCWDGCWWCR
Source	in vitro <i>E.coli</i> expression system
Target Names	AQY2-2
Expression Region	1-149
Protein Length	Full Length
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.