



Recombinant Dog Podoplanin (PDPN)

Product Code	CSB-CF017739DO
Uniprot No.	Q95152
Product Type	Transmembrane Protein
Immunogen Species	Canis lupus familiaris (Dog) (Canis familiaris)
Sequence	ASTVVRPDDIIPGVEDSVVTPGTEDSVVTPGAEDNVVTDGATEEPYESGLTPLVT KNTESVTDLHLEDGPTQESTVHAKEESQSTTTTLNVVTSHSREKVGEDTETTVE KDGLATVTLVGIIVGVLLAIGFIGGIIIVVARKMSGRYSP
Source	in vitro E.coli expression system
Target Names	PDPN
Expression Region	23-169
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.