



Recombinant Human Podoplanin (PDPN)

Product Code	CSB-CF017739HU(A4)
Uniprot No.	Q86YL7
Product Type	Transmembrane Protein
Immunogen Species	Homo sapiens (Human)
Sequence	ASTGQPEDDTETTGLEGGVAMPGAEDDVVTPGTSEDYKSGLTTLVATSVNS VTGIRIEDLPTSESTVHAQEQSPSATASNVATSHSTEKVDGDTQTTVEKDGLST VTLVGIIVGVLLAIGFIGAIIVVMRKMSGRYSP
Source	in vitro E.coli expression system
Target Names	PDPN
Expression Region	23-162
Protein Length	Full Length of Mature Protein
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.