



Recombinant Human Sodium/potassium-transporting ATPase subunit alpha-1 (ATP1A1), partial

Product Code	CSB-BP002322HU
Storage	Store at -20°C, for extended storage, conserve at -20°C or -80°C.
Uniprot No.	P05023
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	<p>TVT VCLTLTAKR MARKNCLVKNLEAVETLGSTSTICSDKTGTLTQNRMTVAHM WFDNQIH EADTTENQSGVSFDKTSATWLALSRIAGLCNRAVFQANQENLPILKRAVAGDA SESALLK CIELCCGSVKEMRERYAKIVEIPFNSTNKYQLSIHKNPNTSEPHLLVMKGAPE RILDRC SSILLHGKEQPLDEELKDAFQNAYLELGGGLGERVLGFCHLFLPDEQFPEGFQF DTDDVNF PIDNLCFVGLISMIDPPRAAVPDAVGKCRSAGIKVIMVTGDHPITAKAIAKGVGII SEGN ETVEDIAARLNIPVSQVNPRDAKACVVHGSCLKDMTSEQLDDILKYHTEIVFAR TSPQQK LIIVEGCQRQGAIVAVTGDGVNDSPALKKADIGVAMGIAGSDVSKQAADMILLD DNFASI VTGVEEGRLIFDNL</p>
Source	Baculovirus
Target Names	ATP1A1
Protein Names	<p>Recommended name: Sodium/potassium-transporting ATPase subunit alpha-1 Short name= Na(+)/K(+) ATPase alpha-1 subunit EC= 3.6.3.9 Alternative name(s): Sodium pump subunit alpha-1</p>
Expression Region	339-772aa
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	partial
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.