




ASY3 Antibody

Product Code	CSB-PA583817XA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0WR66
Immunogen	Recombinant Arabidopsis thaliana ASY3 protein (599-793aa)
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana
Target Names	ASY3

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: ASY3 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 26.2 kDa

Observed band size: 30 kDa