



CpNIFS3 Antibody

Product Code	CSB-PA699834XA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q3E6S9
Immunogen	Recombinant Arabidopsis thaliana CpNIFS3 protein
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order (14-16 weeks)
Target Names	CpNIFS3
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.