




PDLP7 Antibody

Product Code	CSB-PA604499XA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q0WPN8
Immunogen	Recombinant Arabidopsis thaliana PDLP7 protein (31-298aa)
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Target Names	PDLP7

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: PDLP7 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 43.7 kDa

Observed band size: 43.7 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.