




PDF2.2 Antibody

Product Code	CSB-PA658581XA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q39182
Immunogen	Recombinant Arabidopsis thaliana PDF2.2 protein (31-77aa)
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Research Area	Others
Target Names	PDF2.2

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: PDF2.2 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 20.6 kDa

Observed band size: 20.6 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.