



# At3g22700 Antibody

|                     |                                                                                  |
|---------------------|----------------------------------------------------------------------------------|
| Product Code        | CSB-PA885445XA01DOA                                                              |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.                    |
| Uniprot No.         | Q9LUJ1                                                                           |
| Immunogen           | Recombinant Arabidopsis thaliana At3g22700 protein                               |
| Raised In           | Rabbit                                                                           |
| Species Reactivity  | Arabidopsis thaliana                                                             |
| Tested Applications | ELISA, WB (ensure identification of antigen)                                     |
| Form                | Liquid                                                                           |
| Conjugate           | Non-conjugated                                                                   |
| Storage Buffer      | Preservative: 0.03% Proclin 300<br>Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified                                                        |
| Isotype             | IgG                                                                              |
| Clonality           | Polyclonal                                                                       |
| Product Type        | Polyclonal Antibody                                                              |
| Immunogen Species   | Arabidopsis thaliana (Mouse-ear cress)                                           |
| Lead Time           | Made-to-order (14-16 weeks)                                                      |
| Target Names        | At3g22700                                                                        |
| Usage               | For Research Use Only. Not for use in diagnostic or therapeutic procedures.      |