




PEN1 Antibody

Product Code	CSB-PA875527XA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	Q9FR95
Immunogen	Recombinant Arabidopsis thaliana PEN1 protein
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Protein A/G
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Target Names	PEN1

Image



Western Blot

Positive WB detected in Recombinant protein

All lanes: PEN1 antibody at 1:2000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 31.3 kDa

Observed band size: 35 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.