



EIF(ISO)4G2 Antibody

Product Code	CSB-PA250558XA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	O82233
Immunogen	Recombinant Arabidopsis thaliana EIF(ISO)4G2 protein
Raised In	Rabbit
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB (ensure identification of antigen)
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Antigen Affinity Purified
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Lead Time	Made-to-order (14-16 weeks)
Target Names	EIF(ISO)4G2
Usage	For Research Use Only. Not for use in diagnostic or therapeutic procedures.