




KNAT3 Antibody

Product Code	CSB-PA342525XA01DOA
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Uniprot No.	P48000
Immunogen	Recombinant Arabidopsis thaliana KNAT3 protein (108-154aa)
Species Reactivity	Arabidopsis thaliana
Tested Applications	ELISA, WB
Form	Liquid
Conjugate	Non-conjugated
Storage Buffer	Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4
Purification Method	Affinity-chromatography
Isotype	IgG
Clonality	Polyclonal
Product Type	Polyclonal Antibody
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Target Names	KNAT3

Image



Western Blot

Positive WB detected in: recombinant protein

All lanes: KNAT3 Antibody at 1:1000

Secondary

Goat polyclonal to rabbit IgG at 1/50000 dilution

Predicted band size: 23 kDa

Observed band size: 23 kDa

Usage

For Research Use Only. Not for use in diagnostic or therapeutic procedures.