



LTP4 Antibody

| | |
|----------------------------|----------------------------------------------------------------------------------|
| Product Code | CSB-PA873293XA01DOA |
| Storage | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze. |
| Uniprot No. | Q9LLR6 |
| Immunogen | Recombinant Arabidopsis thaliana LTP4 protein |
| Raised In | Rabbit |
| Species Reactivity | Arabidopsis thaliana |
| Tested Applications | ELISA, WB (ensure identification of antigen) |
| Form | Liquid |
| Conjugate | Non-conjugated |
| Storage Buffer | Preservative: 0.03% Proclin 300 Constituents: 50% Glycerol, 0.01M PBS, pH 7.4 |
| Purification Method | Antigen Affinity Purified |
| Isotype | IgG |
| Clonality | Polyclonal |
| Product Type | Polyclonal Antibody |
| Immunogen Species | Arabidopsis thaliana (Mouse-ear cress) |
| Lead Time | Made-to-order (14-16 weeks) |
| Target Names | LTP4 |
| Usage | For Research Use Only. Not for use in diagnostic or therapeutic procedures. |